

PRODUCT INFORMATION AND QUALITY CONTROL SHEET

XYLOSE LYSINE DESOXYCHOLATE (XLD) AGAR

I. INTENDED USE

XLD Agar is a selective and differential medium for the isolation of pathogenic gram negative enteric bacilli, especially *Shigella* species.

II. SUMMARY AND EXPLANATION

Various culture media have been developed for the isolation and differentiation of pathogenic enteric bacilli from sources containing mixed microbiological flora. XLD agar, developed by Taylor, was designed to improve the isolation and identification of enteric pathogens, especially *Shigella*.¹ This formulation improved the growth of more fastidious pathogens that were inhibited by the level of toxic inhibitors incorporated by other selective media. XLD agar differentiates pathogens not only by their nonlactose fermentation but also differentiates them from nonpathogenic, nonlactose fermenters as well.

III. PRINCIPLES OF THE PROCEDURE

The selective agent incorporated in XLD agar, sodium desoxycholate, inhibits the growth of gram positive organisms. Xylose, commonly fermented by most enteric organisms, is included to differentiate *Shigella* species. The shigellae cannot utilize this carbohydrate thus forming red colonies on the medium. *Salmonella* species, which can ferment xylose, would appear indistinguishable from non pathogenic enteric organisms without the presence of lysine. Once the salmonellae deplete the limited supply of xylose in the medium, lysine is utilized via the enzyme lysine decarboxylase. This reaction produces an alkaline by-product that reverses the pH and causes the formation of red colonies. To prevent pH reversion by other non pathogenic enterics, lactose and sucrose were incorporated to produce excess acid.¹

The production of hydrogen sulfide, another differentiating characteristic of the salmonellae, is made evident by the presence of sodium thiosulfate and ferric ammonium citrate. Colonies producing hydrogen sulfide exhibit black centers. Non pathogenic hydrogen sulfide producers, unable to decarboxylate lysine, do not form colonies with black centers because of the excess acid produced by fermenting the alternate carbohydrates lactose and sucrose.

IV. TYPICAL FORMULA AND APPEARANCE

(Approximate formula* per liter of processed water)

Appearance = red, slightly opalescent, no precipitate	
Xylose	3.5 g
L-Lysine	5.0
Lactose	7.5
Sucrose	7.5
Sodium Chloride	5.0
Yeast Extract	3.0
Phenol Red	0.08
Sodium Desoxycholate	2.5
Sodium Thiosulfate	6.8
Ferric Ammonium Citrate	0.8
Agar	13.5

*adjusted and/or supplemented to meet performance criteria.

V. PRECAUTIONS

This product is for IN VITRO diagnostic use only. Culture specimens may contain microorganisms which can be potentially infectious to the user. Strict adherence to aseptic techniques and established precautions against microbiological hazards should be followed throughout the procedure. Carefully dispose of all items which contact patient specimens or isolated bacteria.

VI. STORAGE/SHELF LIFE

Plated media should be stored at 2-8°C (36-46°F), media side

up, in the unopened or resealed package protected from light. DO NOT FREEZE OR EXPOSE TO HIGH TEMPERATURES. Allow unopened plates to warm to room temperature prior to inoculation. Prior to and during inoculation procedures, plates should be handled in a manner that minimizes product exposure to the environment. Product which has exceeded the assigned expiration date noted on the label should not be used.

Do not use plates that exhibit evidence of drying, cracking, discoloration, microbial contamination or any other signs of deterioration. The presence of excessive condensate may indicate plates which have been damaged by exposure to temperature extremes.

VII. SPECIMEN COLLECTION

The quality of culture results depends primarily on the adequacy and condition of the specimen submitted for examination. Proper specimen collection techniques must be followed to ensure the most accurate culture results. Sterile swabs and collection containers should be used. Plates should be inoculated promptly after specimen collection. If a delay in inoculation is unavoidable, transport medium should be employed. Specimens should be collected prior to the initiation of antimicrobial therapy.

Detailed information on proper specimen collection may be obtained from microbiology reference materials.

VIII. MATERIALS PROVIDED

XLD Agar Plates and lot specific Quality Control Certificate.

IX. MATERIALS REQUIRED BUT NOT PROVIDED

Incubator maintaining 33-37°C.
Ancillary culture media, reagents and laboratory equipment as required.

X. PROCEDURE

Inoculate the specimen as soon as possible after it is received in the laboratory. The streak plate method is used primarily to isolate pure cultures from specimens containing mixed flora. If material is being cultured directly from a swab, roll the swab over a small area of the plate surface at the edge (approximately 1/4 to 1/3 of the plate); then streak in a zig-zag fashion with a sterile loop from this inoculated area to cover the entire agar surface. Avoid applying excessive pressure to the agar surface during inoculation to prevent gouging and splitting of the agar media. (Note: Agar surfaces should be smooth and moist but free of excessive moisture which could cause confluent growth patterns.) Incubate plates media side up at 33-37°C for 24-48 hours in an aerobic atmosphere. A nonselective medium should also be inoculated to optimize recovery of enteric pathogenic organisms. This is especially valuable when the total gram negative population is low and to provide an indication of other organisms present in the specimen.

XI. EXPECTED RESULTS

NCCLS Control Organisms (ATCC Strains)

<i>Salmonella typhimurium</i> (ATCC 14028)	Growth, colonies red with black centers
<i>Shigella flexneri</i> (ATCC 12022)	Growth, colonies red
<i>Enterococcus faecalis</i> (ATCC 29212)	Inhibition (partial)
<i>Escherichia coli</i> (ATCC 25922)	Inhibition (partial to complete); colonies yellow to yellow-red

XII. LABORATORY RESULTS

This medium is intended to be used as a primary isolation medium. Presumptive identification of organisms may be made on the basis of typical organism morphology and Gram stain. Definitive identification of organisms and antimicrobial sensitivity determination requires further testing. Additional biochemical information may be obtained from reference microbiology texts.^{2,3,4}

XIII. LIMITATIONS

The ability to detect microorganisms by culture techniques can be affected by the following factors: improper specimen collection, storage and inoculation, initiation of antiinfective therapy prior to specimen collection, improper culture incubation temperatures and atmospheres, improper length of culture incubation, and improper storage and handling of culture media.

XIV. REFERENCES

1. Taylor, W.I. 1965. Isolation of shigellae. I. Xylose lysine agars; new media for isolation of enteric pathogens. *Am. J. Clin. Pathol.*, 44:471-475
2. Finegold, S.M. and W.S. Martin. 1982. *Bailey and Scott's Diagnostic Microbiology*, 6th ed. C.V. Mosby Company, St. Louis.
3. Koneman, E.S., S.D. Allen, V.R. Dowell, Jr. and H. M. Sommers. 1983. *Color Atlas and Textbook of Microbiology*, 2nd ed. J.B. Lippincott Company, Philadelphia.
4. Lennette, E.H., ed. 1985. *Manual of Clinical Microbiology*, 4th ed. American Society for Microbiology, Washington, D.C.

USER QUALITY ASSURANCE/ QUALITY CONTROL PROCEDURES AND INFORMATION

HealthLink recommends that the following quality assurance and quality control procedures be performed on each batch of product.

I. QUALITY ASSURANCE

The following quality assurance procedures must be performed to assure the product will perform according to its intended use within the assigned expiry date:

1. Daily, document that product storage refrigerator maintains temperature within the recommended range: 2-8°C.
2. Daily, document that laboratory incubator maintains temperature within the recommended range: 22-35°C.

II. QUALITY CONTROL

The following incoming inspection procedures must be performed for each batch (batch = same lot, same shipment) of culture media received in the laboratory :

1. Inspect plates according to instructions contained on the "Quality Control Log Sheet." (See also Section VI "STORAGE/SHELF LIFE")
2. Peel off the lower portion of a product bag label (Quality Control Certificate) for the lot being accepted into the laboratory and affix it to the Quality Control Log Sheet.
3. Initial and date the Quality Control Log Sheet.

Note: Notify Technical Service immediately if media does not meet the inspection criteria.

TECHNICAL SERVICE

HealthLink provides a toll free technical service line (1-800-638-2625) to assist with product usage. To receive additional product information, procedural instructions, QA/QC log sheets, Material Safety Data Sheets, or to have technical questions answered, please call between the hours of 9:00 am to 5:00 pm EST.

HealthLink
7779 Bayberry Road
Jacksonville, FL 32256

1-800-638-2625

January, 1997

Product No. 1187 Rev. No. 01