

Biopsy Punches

BIOPUNCH®



SEAMLESS SHARP BLADE

- Stainless steel circular seamless ground blade
- Available in 9 sizes: 1.5, 2, 2.5, 3, 3.5, 4, 5, 6 and 8 mm diameter and a Mix Pack (BP99)










COLOR CODED HANDLES BY SIZE

- BIOPUNCH® handle embossed with size

SAFETY CAPPED FOR YOUR PROTECTION

- Single handed ejection cap
- Clear cap for protection and visual inspection of blade

Ordering information

 # BP15	 # BP40
 # BP20	 # BP50
 # BP25	 # BP60
 # BP30	 # BP80
 # BP35	 # BP99 Mix Pack

Mix Pack contains 5 each of 2, 3, 4, 5 & 6 mm

SURE GRIP HANDLE

- Surgical grip handle design

COLOR CODED DISPENSER BOX

- 25 individual pouches per box
- Reorder level window
- Assorted size box available upon request (Mix Pack contains 5 ea. 2, 3, 4, 5 & 6 mm)

STERILE UNTIL SEAL IS BROKEN

- Temper evident seal - clear seal turns opaque when peeled



For more information, contact your HealthLink sales representative or call **800-638-2625**.
email: HLKCustomerSvc@Clorox.com
visit us: www.hlk.cc
3611 St. Johns Bluff Rd South Ste 1, Jacksonville, FL 32224

© 2014 Clorox Professional Products Company,
1221 Broadway, Oakland, CA 94612

Always use products as directed