

# Ear Curettes



**Ear Curette Child**  
Cat. # 6690 PINK  
2mm



**Ear Curette Spoon**  
Cat. # 6620 GRAY  
2mm



**Ear Curette Oval**  
Cat. # 6610 WHITE  
4mm



**Ear Curette Round**  
Cat. # 6600 TURQUOISE  
3mm

## SAFE, EASY TO USE EAR CURETTE

Our line of ear curettes provides a safe and effective way of removing cerumen (ear wax).

### SAFER THAN STEEL

Steel is very rigid, thus can cause trauma to the ear

### SIMPLER THAN IRRIGATION

Less messy than flushing of ear with water

### REDUCE CHANCE OF CROSS CONTAMINATION

Steel curettes are reusable, thus can create risk of cross contamination

### COLOR CODED

Color coded for easy size and reorder identification  
50 individual curettes per box

